

# **Monitoring of Selected Herbicides, Antibiotics, Steroids, and Industrial Chemicals in Water by ELISA**

**Fernando Rubio, Abraxis LLC**

**Kristy Ramsey, Abraxis LLC**

**Paul Stackelberg, USGS**

**Michael Meyer, USGS**

# Organic Water Contaminants

- Greater Than 90,000 compounds in commerce
- Some manufactured and used in large quantities
- *Examples:* agrochemicals, pharmaceuticals, surfactants, flame retardants, personal care products, commercial and industrial products,
- Usage has increased over the years
- Enter water sources by a variety of environmental pathways

# TARGETED COMPOUNDS

- **Herbicides**  
Acetochlor  
Atrazine
- **Steroid Hormones**  
Estradiol (E2)  
Ethenyl Estradiol (EE2)  
Estrone (E1)
- **Antibiotics**  
Sulfamethazine
- **Surfactants**  
Alkylphenols  
Alkylphenol ethoxylates  
Alkyl ethoxylates
- **Industrial Chemicals**  
Bisphenol A

# Description of Water Samples

## DWT Facility-NJ

- Urban drainage basin serving ~ 850,000 people
- >50 STPs discharge effluent to 2 streams from which DWT withdraws raw water
- Conventional treatment:
  - 1) Screening
  - 2) Addition of powdered AC
  - 3) Addition of acid/base, pH
  - 4) Coagulation
  - 5) Flocculation
  - 6) Sedimentation
  - 7) Primary Desinfection
  - 8) Filtration
  - 9) Secondary Desinfection
  - 10) Add caustic soda, pH

## • Sampling Methods

Twelve water samples collected


Over a 3 week period (July-Aug '03)

Six locations in the facility

- 1) Raw
- 2) Raw and recycled
- 3) Settled
- 4) Disinfected
- 5) Filtered
- 6) Finished

# DWT Facility-NJ

B.



# **Waste Treatment Plant-FL**

## **Purpose of the Project**

**To determine variation of pharmaceuticals, antibiotics, household chemicals, hormones, and industrial waste products going in and coming out of the plan during low and high flow in the wet and dry season.**

# ANALYTICAL METHODS

- Very sensitive methods required sensitivities down to ng/ml, pg/mL)
- HPLC
- LC-MS
- GC-MS
- Immunoassay (ELISA)

# **Immunoassay**

**An immunoassay is an analytical method which uses antibodies as reagents to quantitate specific analytes**




# Immunoassay

- \$6 Billion Industry Worldwide
- 2.5 Billion Tests Sold Annually
- Highly Quantitative
- Regulatory Approved
- Flexible Test Formats
- Diverse Markets and Applications

# Immunoassay Markets


- Clinical
- Agricultural
- Environmental
- Food
- Industrial
- Pharmaceutical
- Veterinary
- Water Quality

# Competitive ELISA




# Competitive ELISA


I. No analyte - high detection signal













II. Analyte present - detection signal reduced



# Competitive ELISA Curve



## ATRAZINE CONCISE FLOWCHART


<p>1.</p>  <p>Separate the rack. Add 200 or 250 <math>\mu</math>L of either Standards, Control or Samples to the bottom of each test tube. <b>NOTE:</b> Same chosen volume should be added for standards and samples.</p>	<p>6.</p>  <p>Invert the combined rack. Blot <b>gently</b>.</p>
<p>2.</p>  <p>Add 250 <math>\mu</math>L of Atrazine Enzyme Conjugate to each test tube.</p>	<p>7.</p>  <p>Add 1 mL of Washing Solution. Wait 2 minutes. Invert the combined rack. Blot <b>gently</b>. Repeat this step.</p>
<p>3.</p>  <p>Add 500 <math>\mu</math>L of mixed Magnetic Particles to each test tube. Vortex.</p>	<p>8.</p>  <p>Separate the rack. Add 500 <math>\mu</math>L of Color Reagent to each test tube. Vortex.</p>
<p>4.</p>  <p>Incubate for 15 minutes.</p>	<p>9.</p>  <p>Incubate for 20 minutes. Prepare blank.</p>
<p>5.</p>  <p>Combine the rack and magnetic base. Seat all tubes. Wait 2 minutes.</p>	<p>10.</p>  <p>Add 500 <math>\mu</math>L of Stopping Solution to each test tube. Read OD 450</p>

For Ordering or Technical Assistance Contact:  
ABRAXIS, LLC 54 Seamwhistle Drive, Warminster, PA 18974  
Phone: 215-357-3911 Fax: 215-357-5232

Atrazine Magnetic Particle Kit Part # 500001, 100 Test





# Microtiter Plate



# SPE Procedure, Hormones




# Herbicide Analysis by ELISA, NJ Samples






# Sulfamethazine Analysis by ELISA, NJ Samples



# E1, E2, EE2 Analysis by ELISA, FL Samples



# E2 ELISA vs. GC, CO Samples

Sample ID	E2, ELISA, (ppt)	E2, GC (ppt)	EE2, ELISA (ppt)	E1, ELISA (ppt)
625A	<0.5	<0.8	1.2	1.4
626A	1.8	2.9	3.9	>50
627A	2.5	2.1	2.8	2.3
628A	<0.5	1.2	5.3	47
629A	13.8	1.4	2.7	18.3
630A	1.8	1.2	2.1	28

# Estradiol Microtiter Plate ELISA



# Atrazine Magnetic Particle Reagents



# Conclusion

- ELISA Technologies, can be used to screen large number of samples timely and cost-effectively.
- Can be use for screening and for quantitation of analytes in water monitoring programs and to efficiently monitor the various stages of treatment at DWT facilities.
- Additional work is currently being performed with surfactants, bisphenol A, as well as comparison to GC values.

# Acknowledgements

- Dr. Larry Barber, USGS Colorado
  - Dr. James Gray, USGS Colorado
- 
- Disclaimer: Mention of a particular product during this presentation does not constitute an endorsement by the USGS.