Department of Defense Environmental Data Quality Issues: Is It Really Perchlorate?





Bill Ingersoll, NAVSEA 04XQ (Labs)
Navy Quality and Accreditation Office
May 2004









Topics



- Data Quality Issues Affecting Missions and Operations
- Perchlorate Testing
 - Current Status of Methods
 - Is It Really Perchlorate?
 - Emerging Technology
- Promoting Sound Science

Perchlorate













Natural and Anthropogenic Sources

- Chilean fertilizer deposits
- New Mexican potash
- Canadian potash
- Californian hanksiteh
- * Bolivian playa crusts

- * Rocket fuel
- Fireworks
- High explosives
- Flares
- * Herbicides
- Automobile airbags
- Tracer munitions
- Detergents

Data Quality Issues Affect Mission and Operations DoD EDQW Charter, DUSD(ES) Sept. 1996:

- Develop DoD policy on sampling, testing & data quality
- Identify Best Practices; standardize across DoD
- Coordinate responses to legislative & regulatory actions
- Represent DoD on environmental data quality Issues
- Respond to DoD IG reports









Data Quality Issues Prescriptive vs. Performance-Based Methods

Office of Water

- Prescriptive methods
- Promulgated (e.g. 40 CFR 136)
- Office of Solid Waster and Emergency Response
 - Performance Based Measurement System (PBMS)
 - SW-846 Methods Manual encourages method modification to meet project needs
 - Performance vs. test procedure oriented
 - Criteria describe data quality needs

State Modifications to PBMS

- Some states legislate SW-846 Methods as requirements
- Performance-based method modifications are not allowed

Data Quality Issues Example: Perchlorate

- Unregulated Contaminant Monitoring Rule (UCMR) List 1 constituent
- Health effects controversial
- Currently 10 states (AZ, CA, FL, IL, MA, MD, NV, NM, NY, and TX) have health advisory limits (1-31 ppb)
- Increasing regulatory pressure to analyze for trace levels (low ppb, sub-ppb)
- Inadequate consideration for complex matrices, method limitations, and measurement uncertainty

Data Quality Issues DoD Interim Policy on Perchlorate, Sept 2003

Components must collect and compile perchlorate data for:

- Drinking water supplies subject to Unregulated Contaminant Monitoring Rule (UCMR)
- Wastewater discharges permitted under National Pollutant Discharge Elimination System (NPDES)
- Environmental restoration sites where release may have occurred and complete human exposure pathway is likely to exist (IR)
- Range assessments conducted under Munitions
 Action Plan and Defense Planning Guidance


Perchlorate Testing Current Status of Methods

- ❖ EPA Method 314.0
 - Ion Chromatography/conductivity detection
 - Drinking Water method
- SW-846 Method 9058 (Nov 2000 Draft)
 - IC/Conductivity detection (with suppressor)
- ❖ PBMS
 - Ion or Liquid Chromatography/Mass Spectrometry (IC or LC/MS)
 - IC or LC/MS/MS
 - Proprietary methods (Modified Method 8321a)


No current, standardized method using MS

EPA Method 314.0 Ion Chromatography/Conductivity

- Currently the only EPAapproved method
- Method Reporting Level(MRL) = 4 ppb(deionized water)
- Compound identification based on retention time
- Conductivity detector is non-specific
 - Responds to any substance that will conduct electric current




EPA Method 314.0 Sensitivity decreases with increasing ionic strength



Is it Really Perchlorate? Example #1

- Treated rocket-motor wastewater, diluted 20X due to high conductivity
- Method 314.0 results:670 ppb(Same result for 2 labs)
- LC/MS/MS results: non-detect






Is it Really Perchlorate Example #1 (cont'd)

To check Method 314:


- 200 ppb perchlorate added to sample
- Spiked peak eluted with 0.6 minute retention-time shift
- Original peak due to interference, not perchlorate
- False positive occurred in every sample from site






Is it Really Perchlorate? Example #2

- Cooling tower water sample diluted 5X because of high conductivity
- Method 314.0 result:194 ppb



Is it Really Perchlorate? Example #2 (cont'd)

- Undiluted sample analyzed by LC/MS/MS
- All four ions detected in proper 3:1 ratio
- LC/MS/MS result:0.98 ppb




Is it Really Perchlorate? Example #3

- Fluidized-bed reactor effluent from perchlorate remediation
- Matrix contains high organic and chloride concentrations
- Method 314.0 result: 220 ppb
- LC/MS/MS result:
- Regulator now willing to accept LC/MS/MS


Is it Really Perchlorate? Example #4

Groundwater Sample

Method 314.0 result: 15 ppb



IC/MS/MS result: 5.64 ppb



Emerging Technology October 2003 Perchlorate Testing Roundtable

- Jointly sponsored by EDQW and IDQTF
- Held in conjunction with EPA Region 6 QA Conference (23 Oct 2003, Dallas)
- Convened technical experts (chemists) from:
 - DoD, EPA and other federal agencies
 - States
 - Academia
 - Analytical laboratories, instrument manufacturers

Emerging Technology Roundtable Objectives

- Discuss current and developing technologies
- Identify most promising technologies
- Recommend path forward for method development

Presentations and Summary can be viewed and downloaded www.navylabs.navy.mil


Emerging Technology Overview

- 'Improved' Ion Chromatography (IC)
- Ion or Liquid Chromatography (IC or LC)
 with single or dual Mass Spectrometry (MS)
- Others
 - Ion-Selective Electrode (ISE)
 - Colorimetry
 - Capillary Electrophoresis
 - Ion Mobility/MS (FAIMS/MS)

Emerging Technology Improved Ion Chromatography

- Preconcentration to improve sensitivity
- Pre-column sample cleanup to reduce ionic strength
- Dual column separation to reduce potential for false positives


Emerging Technology Improved Ion Chromatography (cont'd)





Emerging Technology IC or LC with single or dual Mass Spectrometer

- IC separation & electrospray ionization
- First quadropole
 - MS monitoring99 & 101 m/z
- Collision cell (loses an oxygen)
- Second quadropole
 - MS monitoring83 & 85 m/z




Single MS monitors two ions and one mass ratio

Dual MS monitors four ions and two mass ratios

(Check for 3:1 ion ratio)

LC/MS/MS Sensitivity 0.1 ppb standard in DI water (m/z 83) calculated Method Detection Limit = 0.06ppb)



Performance Summary IC or LC with Single or Dual MS

- Greatly improved perchlorate selectivity and sensitivity
 - Positive identification of perchlorate ion
 - Higher confidence in low-level measurements
 - Sensitivity affected little by changes in ionic-strength (e.g., real world samples)
 - Water taken from the Great Salt Lake and diluted ten fold. Conductivity = 21000 μS MRL 2 ppb.
 - MRL of 0.2 ppb or lower
- Applicable to analysis of complex environmental matrices
- Commercially available as modified SW-846 Method 8321a

Method Performance Comparison Groundwater Samples

Performance Characteristics	Improved IC	IC or LC/MS	IC or LC/MS/MS
Approximate Working range	1-40 ppb	0.1-10 ppb	0.1-10 ppb
MRL	0.5 – 5 ppb*	0.2 ppb	0.2 ppb
Approximate Cost	\$75 - \$125	\$125 - \$150	\$200 - \$250

^{*}dependent on matrix and type of improvement

Promoting Sound Science Roundtable Summary

- Proposed Method 9058 should not be finalized as written
- Method 314.0 modifications can improve sensitivity but not selectivity
- Positive results for perchlorate obtained using Method 314.0 should be confirmed using definitive technology (e.g. MS)

Promoting Sound Science Roundtable Next Steps

- Provide comments on proposed SW-846
 Method 9058 (complete)
- Facilitate process for development/approval of methods employing MS
- Characterize background occurrence of perchlorate
- Pursue development of DoD Guidance for sampling and testing of perchlorate
- Institutionalize the joint roundtable approach for technical problem-solving

Promoting Sound Science EDQW Interim Guidance for Perchlorate Testing

Laboratory requirements:

- Quality systems conformance to DoD QSM
- Nationally-accepted accreditation
- Approval by one or more DoD Component(s)
- Matrix-specific proficiency-testing demonstration

Promoting Sound Science EDQW Interim Guidance for Perchlorate Testing

Method requirements:

- Method Reporting Limit demonstration through analysis of a standard <u>at the MRL</u>
- MRL verification with each batch of samples
- Calibration and quality control that meet or exceed specifications of Method 314.0
- Confirmation of non-specific results (i.e. lon Chromatography) by definitive test methods

Promoting Sound Science Perchlorate Testing Update

DoD:

- EDQW Interim Sampling & Testing Guidance issued by DASN(E) Feb '04
- Perchlorate QAPP Handbook under development
- EPA/ORD & OGW/DW Methods to be finalized Summer '04
- ❖ 314.1 IC/ConductivityMRL 0.5 ug/L

The Payoff

Cost-effective data collection Reproducible data Sound environmental decisions

Improving Environmental Data Quality

... Because the Right Decisions Require Quality Data







SampleJH@navsea.navy.mil www.navylabs.navy.mil www.denix.osd.mil/denix/DOD/Working/EDQW/edqw.html